

# chemagic QA Software

The chemagic QA Software is the ideal solution for labs who require high quality assurance while they use the chemagic Magnetic Separation Module I for nucleic acid isolation. The connection to a standard LIMS system is unproblematic because the QA software generate LIMS compatible log files.

## Features at a glance

- Barcode reading/Sample tracking
- LIMS compatible log files
- Up to date user management
- Easy navigation and set up
- Pre installed protocols
- Tracking the progress of a running protocol



# chemagic Dispenser

The use of the **chemagic MSM I** simplifies and fastens time consuming sample preparation steps like nucleic acid isolation significantly. By using the **chemagic Dispenser** it is even faster and more comfortable. The **chemagic Dispenser** was developed to notably reduce hands-on time preparing a **chemagic MSM I** run. But independent of the amount of needed buffer this particular add-on prefills fast and accurately the tubes/plates during a **chemagic MSM I** run. Variants of the chemagic Dispenser are available fitting to the corresponding chemagic Rod Head formats (12, 24 and 96).

## Features at a glance

- Reduced hands on time
- Simultaneous filling of up to 8 different buffers with different buffer volumes
- Waste of buffer is minimized
- Based on membrane pump technology
- Fully integrated in the chemagic MSM I software
- Easy upgrade of a stand-alone chemagic MSM I (plug and play)



The precise technology of the chemagic Dispenser is also integrated in the chemagic Prepito.